

BTK

Reactivity:Human Mouse

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW:76kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human BTK

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A0500

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:695

Isotype:IgG

Swiss Prot:Q06187

Purity:Affinity purification

For research use only.

To place an order, please [Click HERE](#).