

## RNH1

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:1000 - 1:2000

**Calculated MW:** 50kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human RNH1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

a

**Synonym:**

RAI; RNH;

**Catalog #:** A4079

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 6050

**Isotype:** IgG

**Swiss Prot:** P13489

**Purity:** Affinity purification

For research use only.

**Background:**

Placental ribonuclease inhibitor (PRI) is a member of a family of proteinaceous cytoplasmic RNase inhibitors that occur in many tissues and bind to both intracellular and extracellular RNases (summarized by Lee et al., 1988)

*To place an order, please [Click HERE](#).*