

## EFTUD2

---

**Reactivity:** Human

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 116kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A synthetic peptide of human EFTUD2

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

MFDM; MFDGA; Snu114; Snrp116; SNRNP116; U5-116KD;

**Catalog #:** A4291

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 9343

**Isotype:** IgG

**Swiss Prot:** Q15029

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a GTPase which is a component of the spliceosome complex which processes precursor mRNAs to produce mature mRNAs. Mutations in this gene are associated with mandibulofacial dysostosis with microcephaly. Multiple transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*