

NME2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 17kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NME2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PUF; NDKB; NDPKB; NM23B; NDPK-B; NM23-H2;

Catalog #: A7443

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4831

Isotype: IgG

Swiss Prot: P22392

Purity: Affinity purification

For research use only.

Background:

Nucleoside diphosphate kinase (NDK) exists as a hexamer composed of 'A' (encoded by NME1) and 'B' (encoded by this gene) isoforms. Multiple alternatively spliced transcript variants have been found for this gene. Read-through transcription from the neighboring upstream gene (NME1) generates naturally-occurring transcripts (NME1-NME2) that encode a fusion protein comprised of sequence sharing identity with each individual gene product.

To place an order, please [Click HERE](#).