

LGR5

Reactivity: Human Mouse

Tested applications: WB IHC IF FC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:50 - 1:100 IF 1:20 - 1:100 FC 1:10 - 1:50

Calculated MW: 100kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human LGR5

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FEX; HG38; GPR49; GPR67; GRP49;

To place an order, please [Click HERE](#).

Catalog #: A8090

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8549

Isotype: IgG

Swiss Prot: O75473

Purity: Affinity purification

For research use only.