

NKX3-2

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human NKX3-2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NK3 homeobox 2; BAPX1, Homeobox protein Nkx 3.2, NK3 homeobox 2, NKX3 2, NKX3.2, NKX3-2, NKX3B;

Background:

This gene encodes a member of the NK family of homeobox-containing proteins. The encoded protein may play a role in skeletal development.

To place an order, please [Click HERE](#).

Catalog #: A8093

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 579

Isotype: IgG

Swiss Prot: P78367

Purity: Affinity purification

For research use only.