

CNTROB

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:101kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CNTROB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LIP8; PP1221;

Catalog #:A8277

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:116840

Isotype:IgG

Swiss Prot:Q8N137

Purity:Affinity purification

For research use only.

Background:

This gene encodes a centrosomal protein that interacts with BRCA2, and is required for centriole duplication and cytokinesis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).