

GTF2H2C_2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 44kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GTF2H2C_2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GTF2H2D;

To place an order, please [Click HERE](#).

Catalog #: A8297

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 730394

Isotype: IgG

Swiss Prot: Q6P1K8

Purity: Affinity purification

For research use only.