

SAE1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:500

Calculated MW: 38kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SAE1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9960

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10055

Isotype: IgG

Swiss Prot: Q9UBE0

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).