

CEP164

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 164kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CEP164

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A9964

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22897

Isotype: IgG

Swiss Prot: Q9UPV0

Purity: Affinity purification

For research use only.