

alpha Tubulin

Reactivity: Human Mouse

Tested applications: WB

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human alpha Tubulin

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: AC025

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10376

Swiss Prot: P68363

For research use only.