

GST-Tag

Tested applications:WB IF IP

Recommended Dilution:WB 1:2000 - 1:5000 IF 1:50 - 1:200 IP 1:20 - 1:50

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of GST tag

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:AE001

Antibody Type:

Monoclonal Antibody

Species:Mouse

Isotype:IgG

Purity:Affinity purification

For research use only.