

TAP-tag

Tested applications:WB IP

Recommended Dilution:WB 1:1000 - 1:5000 IP 1:100 - 1:500

Calculated MW:48kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of TAP epitope

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:AE015

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Isotype:IgG

Purity:Affinity purification

For research use only.

To place an order, please [Click HERE](#).