

Phospho-Lck-Y416/FYN-Y416

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 58kDa

Observed MW: Refer to figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y416 of human Lck and Y416 of human FYN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bd

Catalog #: AP0491

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3932; 2534

Isotype: IgG

Swiss Prot: P06239; P06241

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).