

IFN a 1b Human

Description: Interferon-alpha 1b Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 166 amino acids and having a molecular mass of 19392 Dalton. The Interferon-alpha 1b gene was obtained from human leukocytes. The IFN-a 1b is purified by proprietary chromatographic techniques.

Catalog #: CYPs-290

For research use only.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Amino Acid Sequence: The sequence of the first five N-terminal amino acids was determined and was found to be Cys-Asp-Leu-Pro-Glu.

Purity: Greater than 96.0% as determined by: (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

Lyophilized from a 0.2

Stability:

Lyophilized Interferon alpha 1b although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution IFN-alpha 1b should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized IFN alpha 1b in sterile 18M-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Introduction:

IFN-alpha 1b is produced by macrophages, IFN-alpha 1b has antiviral activities and stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

Biological Activity:

The specific activity as determined in a viral resistance assay using bovine kidney MDBK cells was found to be 1 x 10⁸ IU/ mg.

To place an order, please [Click HERE](#).