

Enterokinase Porcine

Description: a fusion tag is located in the N-terminus with an enterokinase site, enterokinase will be able to remove the fusion tag and to generate the protein exactly as you need without adding any unwanted residues. NeoBiolabs enterokinase is a highly purified enterokinase from porcine.

The enzyme has been extensively purified and tested to ensure that there are no other contaminating proteases.

Catalog #:ENPS-274

For research use only.

Synonyms:Enteropeptidase, EC 3.4.21.9, Enterokinase, Serine protease 7, ENTK, MGC133046.

Source:Porcine.

Physical Appearance:Sterile Liquid.

Formulation:

2 IU/

Stability:

One year when stored at -20°C, one week at room temperature.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Enteropeptidase or enterokinase is an enzyme involved in human digestion. It is produced by cells in the duodenum wall, and is secreted from duodenum's glands, the crypts of Lieberkhn, whenever ingested food enters the duodenum from the stomach. Enteropeptidase has the critical job of turning trypsinogen(a zymogen) to trypsin, indirectly activating a number of pancreatic digestive enzymes. Enteropeptidase is a serine protease enzyme (EC3.4.21.9). Enteropeptidase is a part of the Chymotrypsin-clan of serine proteases, and is structurally similar to these proteins.

To place an order, please [Click HERE](#).