

ABHD10 Human

Description: ABHD10 Human Recombinant produced in E. coli is a single polypeptide chain containing 279 amino acids (53-306) and having a molecular mass of 30.9kDa. ABHD10 is fused to a 25 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: ENPS-619

For research use only.

Synonyms: Abhydrolase domain containing 10 mitochondrial2, FLJ11342, EC 3.4.-.-.

Source: E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHH SSGLVPRGSH MGSMTLSLF LNRDLPNLA
YKKLKGKSPG IIFIPGYLSY MNGTKALAE EFCKSLGHAC IRFDYSGVGS SDGNSEESTL
GKWRKDVLSI IDDLADGPQI LVGSSLGGWL MLHAAIARPE KVVALIGVAT AADTLVTKFN
QLPVELKKEV EMKGVWSMPS KYSEEGVYNV QYSFIKEAEH HCLLHSPIPV NCPIRLHGM
KDDIVPWHTS MQ

Purity: Greater than 90% as determined by SDS-PAGE.

Formulation:

The ABHD10 solution (1mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 100mM NaCl and 20% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Abhydrolase domain-containing protein 10 (ABHD10) is a member of the AB hydrolase superfamily.

To place an order, please [Click HERE](#).