

H1N1 Solomon Islands

Description: Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza A virus, strain A/Solomon Islands/03/06. The Influenza Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient.

Catalog #: IHPS-029

For research use only.

Physical Appearance: Sterile Filtered colorless solution.

Purity: Greater than 90.0% as determined by SDS-PAGE.

Formulation:

The H1N1 A/Solomon Islands/03/06 solution contains 0.1M NaCl, 10mM Tris Hcl, 1mM EDTA pH-8(STE), 0.1 % sodium azide (NaN₃) and 0.005 % thimerosal.

Stability:

A/Solomon Islands/03/06 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

H1N1 is a subtype specie of Influenza A virus. H1N1 Influenza Virus has mutated into various strains such as the Spanish Flu strain, mild human flu strains, endemic pig strains, and various strains found in birds. The Influenza A Virus is a globular particle about 100nm in diameter, sheathed in a lipid bilayer derived from the plasma membrane of its host. Studded in the lipid bilayer are two integral membrane proteins some 500 molecules of hemagglutinin ("H") and some 100 molecules of neuraminidase ("N"). Within the lipid bilayer are 3000 molecules of matrix protein and 8 pieces of RNA. Each of the 8 RNA molecules is associated with many copies of a nucleoprotein, several molecules of the three subunits of its RNA polymerase some "non-structural" protein molecules of uncertain function.

To place an order, please [Click HERE](#).