

ANXA7 Human

Description:ANXA7 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 490 amino acids (1-466 a.a.) and having a molecular mass of 52.9kDa. ANXA7 is fused to a 24 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #:PRPS-456

For research use only.

Synonyms:Annexin A7, Annexin-7, Annexin VII, Synexin, ANXA7, ANX7, SNX.

Source:Escherichia Coli.

Physical Appearance:Sterile Filtered clear solution.

Amino Acid Sequence:MGSSHHHHHH SSGLVPRGSH MGSMSYPGY PPTGYPPFPG
YPPAGQESSF PPSGQYPYPS GFPPMGGGAY PQVPSSGYPG AGGYPAPGGY PAPGGYPGAP
QPGGAPSYPG VPPGQGFGVP PGGAGFSGYP QPPSQSYGGG PAQVPLPGGF
PGGQMPSQYP GGQPTYPSQP ATVTQVTQGT IRPAANFDAI RDAEILRKAM KGFGTDEQAI
VDVVANRSND QRQKIAAFK TS

Purity:Greater than 85.0% as determined by SDS-PAGE.

Formulation:

ANXA7 protein solution (0.25mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 150mM NaCl, 1mM DTT and 40% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Annexin VII is a member of the annexin family of calcium-dependent phospholipid binding proteins. Annexin VII has molecular weight of approximately 51 kDa with a unique, highly hydrophobic N-terminal domain of 167 amino acids and a conserved C-terminal region of 299 amino acids. The latter domain is composed of alternating hydrophobic and hydrophilic segments. Structural analysis of the protein suggests that Annexin VII is a membrane binding protein with diverse properties including voltage-sensitive calcium channel activity, ion selectivity and membrane fusion.

To place an order, please [Click HERE](#).