

CD45 Human

Description: CD45 Human Recombinant (aa 1102-1150) expressed in E.coli, shows a 33 kDa band on SDS-PAGE (Including GST). The CD45 is purified by proprietary chromatographic techniques.

Catalog #: PRPS-489

For research use only.

Synonyms: Leukocyte common antigen, EC 3.1.3.48, L-CA, T200, CD45 antigen, PTPRC, LCA, LY5, B220, CD45, GP180.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Formulation:

CD45 at 100

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

ELISA. Inhibition Assays. Western Blotting.

Introduction:

CD45 leukocyte common antigen (LCA) belongs to the family of at least four isoforms of membrane glycoproteins (220, 205, 190, 180kDa) expressed on hematopoietic cell lines but absent on non-hematopoietic cell lines, normal and malignant non-hematopoietic tissues. The intracellular portion of these molecules has protein phosphatase activity and is involved in regulation of transmembrane signals.

To place an order, please [Click HERE](#).